

# Buffer Compatibility with Olink

Buffers can be used for preparation of cell and tissue lysates, elution of proteins from solid supports, and lavage solutions. The main considerations for buffer formulation are: i) proteins should be maintained under native conditions for proper antibody recognition, and ii) be careful of reagents that can interfere with PCR.

## Recommendations for Buffer Formulation

---

### Detergents

- Excessive detergent concentration can lead to protein denaturation
- The concentration of ionic detergents should be  $\leq 0.1\%$ , including sodium dodecyl sulfate (SDS) and deoxycholic acid (DCA; also known as deoxycholate)
- Non-ionic detergent concentrations, such as Tween® 20, Triton™ X100, and NP-40 should be  $\leq 1\%$

### Other denaturants to avoid

- Dithiothreitol (DTT) in excess of 1 mM
- Alcohol
- Heat
- Urea
- Heavy metal salts

### Salts

- High salt concentrations can promote the aggregation and precipitation of proteins
- Concentrations of common salts should be:  $\leq 250$  mM NaCl,  $\leq 25$  mM KCl, and  $\leq 10$  mM  $\text{MgCl}_2$

### pH

- The pH of buffers should be close to physiological levels (7-8 range)

### Commercial buffers

- T-PER™ Tissue Protein Extraction Reagent (#78510) and M-PER™ Mammalian Protein Extraction Reagent (#78501) from Thermo Fisher are recommended.
- Bio-Plex® Cell Lysis Kit (Bio-Rad #171304011). Add the following to the buffer just prior to use: Factors 1 and 2 (from the kit) + protease inhibitor cocktail.
- RIPA buffer can be custom made as: 50 mM Tris-HCl (pH 7.4), 150 mM NaCl, 1 mM EDTA, 1% Triton X100, 0.1% DCA. We have tested commercial brands, and Millipore RIPA Lysis Buffer (#20-188) works well with Olink technology, whereas Sigma RIPA Buffer (#R0278) does not.
- NP40 lysis buffer can be custom made as: 1X Tris-EDTA buffer (pH 8.0), 1% NP-40, 0.1% Triton X100, 0.1% sulfobetaine, 150 mM NaCl, and 1X protease inhibitor cocktail.

## Diluents

- The original lysis/elution buffer should be used when normalizing samples to a set protein concentration
- Olink Diluent should be used when Olink panels or alternative matrices require pre-dilution of samples

## PCR inhibitors

- EDTA should be  $\leq 25$  mM to avoid sequestration of  $Mg^{2+}$
- Urea (found in urine; pre-dilute samples at least 1:4 with Olink Diluent)
- Heme (found in hemolysates; refer to validation documents for acceptable limits)

## Fluorescent dyes

- When running samples on our Target platform, be careful of the presence of fluorescent dyes that are known to interfere with quantitative PCR, such as fluorescein

## Protease and phosphatase inhibitors

- It is generally recommended to include protease inhibitors within cell and tissue lysis buffers
- Roche cOmplete™ Mini Protease Inhibitor Cocktail (#11836153001) is highly recommended. One tablet can be dissolved in 10 ml of lysis buffer. Alternatively, a 10X solution can be prepared by dissolving 1 tablet in 1 ml of distilled water or PBS, or a 7X stock in 1.5 ml. The stock solution can be stored at 4°C for  $\leq 2$  weeks or -20°C for  $\leq 12$  weeks. Use a 1X final concentration of inhibitor cocktail and avoid excess final concentrations (e.g., 2X or 3X).
- General concentrations of protease inhibitors known to be compatible are: 1 mM PMSF, 10 mM AEBSF, 8 mM aprotinin, 0.2 mM leupeptin, 0.4 mM bestatin, 0.15 mM pepstatin, and 0.15 mM E-64
- Olink assays do not recognize phosphorylation sites, but phosphatase inhibitors such as NaF in the range of 5-10 mM and  $Na_3VO_4$  at 1-2 mM are acceptable, as well as cocktails such as PhosSTOP™

## Multi-omic studies

- Protein samples for multi-omics studies should be processed separately since many protocols for RNA/DNA isolation contain strong denaturants such as phenol and/or guanidinium thiocyanate
- PAXgene® tubes, Allprotect® Tissue Reagent, and RNeasy® are not compatible with Olink
- Cell-Free DNA BCT® tubes from Streck are compatible with Olink

Please contact [support@olink.com](mailto:support@olink.com) for further information on running alternative matrices.

[www.olink.com](http://www.olink.com)

For research use only. Not for use in diagnostic procedures.  
This product includes a license for non-commercial use. Commercial users may require additional licenses. Please contact Olink Proteomics AB for details.  
There are no warranties, expressed or implied, which extend beyond this description.  
Olink Proteomics AB is not liable for property damage, personal injury, or economic loss caused by this product.  
OLINK, NPX, and the Olink logotype are trademarks registered, or pending registration, by Olink Proteomics AB. All third-party trademarks are the property of their respective owners.  
Olink Proteomics, Dag Hammarskjölds väg 52B, SE-752 37 Uppsala, Sweden  
AM-LEB, v1.3